

PRODUCT OVERVIEW

# Aurora<sup>®</sup> Rapid<sup>™</sup> 8x75

8 cm x 75 µm



## Maximal throughput for single-cell and low input samples.

Aurora<sup>®</sup> Rapid<sup>™</sup> 8x75 columns deliver high-throughput analysis of single cells and small sample loads, providing extreme sensitivity in shorter gradients. With further improved spray stability and increased robustness in Generation 4 Aurora<sup>®</sup> Rapid<sup>™</sup> 8x75, achieve exceptional depth of coverage and high throughput without compromising data quality.

### Product Benefits

- + High throughput
- + Nanoflow
- + Sensitivity
- + Analyse more cells in less time

Ideal for:  
High-throughput  
single cell  
proteomics


50-100  
SPD

10-30<sub>min</sub>  
Gradients

### Specifications

<b>Column format</b>	Analytical column
<b>Column type</b>	Reversed-phase
<b>For use with</b>	UHPLC
<b>Length</b>	8 cm
<b>Inner Diameter</b>	75 µm
<b>Pore size</b>	120 Å
<b>Pressure</b>	>1700 bar
<b>Temp. limits</b>	60°C
<b>Particle size</b>	1.7 µm
<b>pH stability</b>	1-8
<b>Stationary phase</b>	C18
<b>Suggested flow rate</b>	100-1000 nL/min


### Compatibility



	Aurora <sup>®</sup> CSI	Aurora <sup>®</sup> XT
CaptiveSpray	✓	
CaptiveSpray 2	✓	
CaptiveSpray Ultra	✓	
EASY-Spray		✓
Nanospray Flex		✓
Newomics UniESI		✓
Newomics DuoESI		✓
OptiFlow Turbo V		
OptiFlow Pro		





Introducing the HeatSync<sup>™</sup> range - the perfect companion to make the world's best chromatography even better. Includes the new HeatSync<sup>™</sup> Controller and HeatSync<sup>™</sup> Column Heater.



**Ion Opticks Pty Ltd**  
 ABN: 99 621 674 459  
 12 Gipps St  
 Collingwood VIC 3066  
 Australia

[www.ionopticks.com](http://www.ionopticks.com)

## More than 3200 proteins identified from inputs as low as 125 pg




**Figure 1:** Unique protein group identifications across different columns. A HeLa tryptic digest (250 pg) was separated on an Aurora Rapid 8 cm × 75 μm column using an 80 SPD method on a Vanquish Neo system coupled to a Bruker timsTOF Ultra 2 mass spectrometer. Three runs were performed per column. Raw data were analyzed in Spectronaut version 19.8 with the “Match between runs” feature disabled

**Figure 2:** Unique peptide identifications across different columns. A HeLa tryptic digest (250 pg) was separated on an Aurora Rapid 8 cm × 75 μm column using an 80 SPD method on a Vanquish Neo system coupled to a Bruker timsTOF Ultra 2 mass spectrometer. Three runs were performed per column. Raw data were analyzed in Spectronaut version 19.8 with the “Match between runs” feature disabled.



## Incredibly stable retention times

**Figure 3:** Stable peptide retention times across different columns. Ten peptides were selected and their retention times assessed across all columns from 3 HeLa tryptic digest injections (250 pg) on Aurora Rapid 8 cm × 75 μm columns, using an 80 SPD method. Samples were run on a Vanquish Neo and Bruker timsTOF Ultra 2 mass spectrometer. Raw data were analyzed in Spectronaut version 19.8 with the “Match between runs” feature disabled.



## Column backpressure stability and protein group IDs over extended use

**Figure 4:** Backpressure stability and protein group IDs over extended column use from HeLa tryptic digest QC injections (250 pg) on Aurora Rapid 8 cm × 75 μm columns, using a 60 SPD method. Samples were run on a Vanquish Neo and Bruker timsTOF Ultra 2 mass spectrometer. Raw data were analyzed in Spectronaut version 19.8 with the “Match between runs” feature disabled.



Further literature at:  
[ionopticks.com](http://ionopticks.com)  
**ionopticks**